



Port list

No.	Size	Connection	Usage
N1	1.5 ^s (ϕ 38)	Tri-clamp	Disinfection water outlet
N2	1.5 ^s (ϕ 38)	Tri-clamp	water inlet
N3	ϕ 350		manhole
N4	1.25 ^s (ϕ 32)		breath port
N5	1.5 ^s (ϕ 38)	Tri-clamp	CIP
N6	1.5 ^s (ϕ 38)	Tri-clamp	return port
N7	1 ^s (ϕ 25)	Tri-clamp	level gauge

Technical requirement

1. Total volume: 630L, effective volume: 500L;
2. Protect the board surface inside and outside, and throw a fine sanding tape at the weld;
3. Welding adopts argon arc welding;
4. After making, fill in water to test for leaks.

借(通)用
件 登 记
旧底图总号
底图总号
签 字
日 期
档案员 日期

